

SOP – Thermal Cyclers (Bio-Rad)

Contacts		
Emergency: Dan Voicu Axel Guenther	dan.voicu@utoronto.ca guenther@mie.utoronto.ca	(416)-978-8473 (416)-978-1282
Training: Victor Sit Dan Voicu	victor.sit@utoronto.ca dan.voicu@utoronto.ca	[TBD] (416)-978-8473


Equipment Overview

The Bio-Rad thermal cyclers are touchscreen devices used for PCR experiments to amplify RNA or DNA sequences. However, other experiments can be programmed as long they can fit into the typical 0.2 mL tubes (T100) or plates (96-/384-well plates). There are 3 types of instruments at the facility:

Instrument	Sample types	Features	Max Ramp Rate
CFX96	96-well plate	Incubation, real-time PCR with 6 filtered LEDs (450-730 nm)	5°C/sec
CFX384	384-well plate	Incubation, real-time PCR with 5 filtered LEDs (450-690 nm)	2.5°C/sec
T100	0.2 uL tubes, 96 well-plate	Incubation, PCR	4°C/sec

Notes before use:

- Ensure plate is clean and dry before placing it into the thermocycler
- Only use 70% ethanol for cleaning to prevent corrosion from chemicals.
- The temperature range 4–100°C.
- Book the instruments ahead of time on LMACs.
- Warning: the heating block and lid can be hot during operation.
- Below is the list of fluorophores available for the CFX96/CFX384:



Supplies and Reagents

Provided by CRAFT:


- 70% Ethanol
- Kimwipes
- PPE (lab coat, gloves)
- Plastic insert (for tubes on the T100) – one side is for tubes with dome caps, the other is for flat caps.




Method

Setup:


1. Log into LMACs to turn on the thermal cycler.
2. On the home screen, select **Saved Files**:




3. Create a folder for your lab (ie. Professor's surname, PI) under **CT003109**:



4. Once a folder has been created, click the **Home** button, and select **New Protocol**.
5. Create a protocol for your assay.




6. Click **Save** and select the designated folder (I.e. Professor's surname, PI):



7. Click **Run** when ready to start the experiment.

8. After the experiment is finished, ensure to cancel or stop the run before shutting down the system.

9. Export the data by selecting **Export** and plugging an external drive to the USB A port:



10. Remove your sample, and log out of LMACs.

Cleanup:

1. Close the software.

2. Wipe the workstation with RTU or 70% ethanol.

Troubleshooting:

1. If the touchscreen is frozen, restart the instrument.
2. If issues cannot be resolved, contact a CRAFT staff member.

Resources:

T100 specifications: https://www.bio-rad.com/webroot/web/pdf/lsr/literature/Bulletin_6060.pdf

CFX96 specifications: https://www.bio-rad.com/webroot/web/pdf/lsr/literature/Bulletin_6093.pdf

CFX384 specifications: https://www.bio-rad.com/webroot/web/pdf/lsr/literature/Bulletin_6096.pdf

User Manual (CFX96/CFX384): <https://www.bio-rad.com/en-ca/node/39981>